

LED Facial Mask + Neck Quick Start

A simple setup guide for using your Alpha LED facial mask and neck piece as part of a calm at-home beauty routine.



Quick Start

- 1 Start with clean, dry skin. Remove creams or products unless your product guide says they are suitable.
- 2 Check the mask, neck piece, cable, controller, and connection points before use.
- 3 Position the mask comfortably and secure the straps gently. Do not pull it tight.
- 4 Use the timing, mode, and brightness settings shown in your product guide. Start gently if you are new to LED devices.
- 5 After use, switch off, disconnect, and wipe the mask with a soft dry cloth before storing it safely.

Care + Safety

- Avoid looking directly into active LEDs. Use eye protection or follow the eye-safety instructions supplied with your model.
- Do not use on broken, irritated, or unusually sensitive skin.
- Keep the device dry and do not submerge it in water.
- Stop if anything feels uncomfortable.

Important This guide is general product support only and is not medical advice. Results and experiences vary.